


THE WRAP SERIES



Single Vision Spherical

Designed For Wrap Frames

SuperLite™ 1.60
Coated and Uncoated



Index Nd	1.592
Abbe Value Nd	42.0
Density gr/cm ³	1.30

Blank Specifications

Nominal Base	True Base 1.530	Base Radius [mm]	Sag 50mm*	Concave Curve	CT [mm]	ET [mm]	Diameter	Power Range
6.00	5.89	89.98	3.54	6.10	10.0	10.4	80	-9.00 to +5.75
7.00	6.65	79.66	4.02	6.10	11.1	10.3	80	-8.50 to +6.50
8.00	7.86	67.43	4.81	8.25	10.0	11.0	80	-6.75 to +8.25

Prescription Range - Sphere: -9.00[D] to +8.25[D] Cylinder: up to -6.00[D]

Minimum center for surfacing 1.5mm