


THE WRAP SERIES

Single Vision Spherical

For Wrap Frames

SuperLite™ 1.67
Transitions® XTRActive™



Index Nd	1.664
Abbe Value Nd	32.0
Density gr/cm ³	1.37

Blank Specifications

Nominal Base	True Base 1.530	Base Radius [mm]	Sag 50mm*	Concave Curve	CT [mm]	ET [mm]	Diameter	Power Range
6.00	5.89	89.98	3.54	6.10	10.0	10.3	77	-10.25 to +6.50
7.00	6.69	79.22	4.05	6.10	11.5	10.7	77	-9.50 to +7.50
8.00	7.86	67.43	4.81	8.25	9.7	10.7	77	-7.50 to +9.25

Prescription Range - Sphere: -10.25[D] to +9.25[D] Cylinder: up to -6.00[D]

Minimum center for surfacing 1.5mm

