


THE WRAP SERIES

Single Vision Spherical

Trivex For Wrap Frames
Coated and Uncoated

Index Nd	1.53
Abbe Value Nd	43.0
Density gr/cm ³	1.11



Blank Specifications

Nominal Base	True Base 1.530	Base Radius [mm]	Sag 50mm*	Concave Curve	CT [mm]	ET [mm]	Diameter	Power Range
6.25	6.38	83.07	3.851	6.25	14.0	14.0	76	-6.00 to +6.00
7.25	7.46	71.05	4.544	6.25	16.1	14.0	76	-6.00 to +6.00
8.25	8.45	65.72	5.198	6.25	17.8	14.0	76	-6.00 to +6.00

Prescription Range - Sphere: -6.00 [D] to +6.00[D] Cylinder: up to -6.00[D]

Minimum center for surfacing 1.5mm