




Shamir SV™ with Attitude

Polycarbonate 1.586

Polarized Gray



Index Nd	1.586
Abbe Value Nd	30.0
Density gr/cm ³	1.20

Blank Specifications

Nominal Base	True Base 1.530	Base Radius [mm]	Sag 50mm*	Concave Curve	CT [mm]	ET [mm]	Diameter	Power Range
8.00	8.29	63.93	5.09	9.50	10.85	13.1	76	-4.00 to +2.00

Prescription Range - Sphere: -4.00[D] to +2.00[D] Cylinder: up to -4.00[D]

Minimum center for surfacing 1.5mm