


SuperLite™ 1.67

Clear

Coated and Uncoated

Material Data:

Index Nd	1.664
Abbe Value Nd	32.0
Density gr/cm ³	1.37



Minimum Fitting Height: 19mm	PRP Decentration: 2.5mm	Fitting Point: 4mm above PRP
------------------------------	-------------------------	------------------------------

*** Do not sag! Use printed sag values for power calculations only.
 Minimum Center Thickness for Surfacing 1.5mm.**

Blank Specifications

Nominal Base	True Base 1.530	Base Radius [mm]	Sag 50mm*	Concave Curve	CT [mm]	ET [mm]	Diameter	Add Powers
1.50	1.51	352.04	0.89	6.10	11.6	17.1	75 / 80	1.00 – 3.50
3.00	3.15	168.07	1.87	6.10	10.2	13.6	75 / 80	1.00 – 3.50
5.00	5.18	102.33	3.10	6.10	7.7	8.2	75 / 80	1.00 – 3.50
7.50	7.31	72.48	4.45	6.10	10.9	8.1	75 / 80	1.00 – 3.50

Prescription Range - Sphere: -12.50[D] to +7.50[D] Cylinder: up to -6.00[D]

Nominal Base	Power Range
1.50	-12.50 to -6.25
3.00	-6.00 to -2.00
5.00	-1.75 to +2.25
7.50	+2.50 to +7.00