


Trivex™ 1.53

Coated and Uncoated

Material Data:

Index Nd	1.533
Abbe Value Nd	45.0
Density gr/cm ³	1.31



Ink Marking: Black. Invisible Marking: Red

Minimum Fitting Height: 19mm	PRP Decentration: 2.5mm	Fitting Point: 4mm above PRP
------------------------------	-------------------------	------------------------------

*** Do not sag! Use printed sag values for power calculations only.
 Minimum Center Thickness for Surfacing 1.5mm.**

Blank Specifications

Nominal Base	True Base 1.530	Base Radius [mm]	Sag 50mm*	Concave Curve	CT [mm]	ET [mm]	Diameter	Add Powers
3.00	3.18	166.67	1.89	4.32	9.0	10.0	76/81	1.00 - 3.00
5.00	5.24	101.15	3.14	6.34	7.0	7.8	76/81	1.00 - 3.00
7.50	7.34	72.21	4.47	6.34	10.8	7.8	76/81	1.00 - 3.00

Prescription Range - Sphere: -5.75[D] to +6.00[D] Cylinder: up to -6.00[D]

Nominal Base	Power Range
3.00	-5.75 to -2.00
5.00	-1.75 to +2.25
7.50	+2.50 to +6.00