




# SuperLite™ 1.67

## Coated and Uncoated

### Material Data:

Index Nd	1.664
Abbe Value Nd	32.0
Density gr/cm <sup>3</sup>	1.37



Ink Marking: Black. Invisible Marking: Red

Minimum Fitting Height: 16mm	PRP Decentration: 2.5mm	Fitting Point: 4mm above PRP
------------------------------	-------------------------	------------------------------

**\* Do not sag! Use printed sag values for power calculations only.  
Minimum Center Thickness for Surfacing 1.5mm**

### Blank Specifications

Nominal Base	True Base 1.530	Base Radius [mm]	Sag 50mm*	Concave Curve	CT [mm]	ET [mm]	Diameter	Add Powers
1.50	1.48	357.27	0.88	6.10	10.2	15.8	72/77	1.00 – 3.00
3.00	2.98	178.07	1.76	6.10	8.8	12.5	72/77	1.00 – 3.00
4.50	4.51	117.42	2.69	6.10	9.1	10.8	72/77	1.00 – 3.00
5.75	5.75	92.23	3.45	6.10	8.2	8.4	72/77	1.00 – 3.00
7.00	7.00	75.66	4.25	6.10	10.1	8.4	72/77	1.00 – 3.00

**Prescription Range - Sphere: -12.50[D] to +6.00[D] Cylinder: up to -6.00[D]**

Nominal Base	Power Range
1.50	-11.25 to -6.25
3.00	-6.00 to -2.00
4.50	-1.75 to +2.00
5.75	+2.25 to +4.00
7.00	+4.25 to +6.00