




SuperLite™ 1.6
Transitions Gray



Material Data:

Index Nd	1.592
Abbe Value Nd	42.0
Density gr/cm ³	1.30

Ink Marking: Black. Invisible Marking: Red

Minimum Fitting Height: 16mm	PRP Decentration: 2.5mm	Fitting Point: 4mm above PRP
------------------------------	-------------------------	------------------------------

*** Do not sag! Use printed sag values for power calculations only.**
Minimum Center Thickness for Surfacing 1.5mm

Blank Specifications

Nominal Base	True Base 1.530	Base Radius [mm]	Sag 50mm*	Concave Curve	CT [mm]	ET [mm]	Diameter	Add Powers
1.50	1.48	356.92	0.88	6.10	9.7	14.6	72/77	1.00 – 3.50
3.00	3.13	169.12	1.86	6.10	10.3	13.1	72/77	1.00 – 3.50
5.00	5.16	102.77	3.09	6.10	8.0	8.1	72/77	1.00 – 3.50
7.50	7.32	72.45	4.45	6.10	11.1	8.1	72/77	1.00 – 3.50

Prescription Range - Sphere: -11.25[D] to +7.00[D] Cylinder: up to -6.00[D]

Nominal Base	Power Range
1.50	-11.25 to -6.25
3.00	-6.00 to -2.00
5.00	-1.75 to +2.25
7.50	+2.50 to +7.00