


Single Vision HC & UC Spherical



1.60 Clear

Certified for Shamir FREEFORM®

Index Nd	1.596
Abbe Value	38
Density gr/cm ³	1.31



NOMINAL BASE CURVE	TRUE BASE (1.530)	BASE RADIUS (MM)	O.SAG (50MM)	CONCAVE CURVE	CT (MM)	ET (MM)	DIAMETER
0.50	0.51	1036.21	0.30	5.87	8.20	15.90	76
2.00	1.87	283.49	1.10	7.73	8.00	16.90	76
2.75	2.80	188.99	1.66	7.73	6.40	13.90	76
3.75	3.75	141.41	2.23	7.73	7.80	14.00	76
5.00	5.04	105.25	3.01	5.87	8.20	9.60	76
6.00	6.03	87.87	3.63	5.87	9.00	9.00	76
7.00	7.02	75.50	4.25	5.87	10.20	8.40	76
8.00	8.02	66.12	4.91	5.87	10.00	6.60	76
9.00	9.01	58.81	5.57	7.73	13.20	11.20	76

Minimum Center Thickness for Surfacing 2.0mm



SHAMIR

shamirlens.com • 877.514.8330