


Single Vision HC UC Spherical



1.67 Clear

Certified for Shamir FREEFORM®

Index Nd	1.66
Abbe Value	32
Density gr/cm ³	1.36



NOMINAL BASE CURVE	TRUE BASE (1.530)	BASE RADIUS (MM)	SAG (50MM)	CONCAVE CURVE	CT (MM)	ET (MM)	DIAMETER
0.50	0.50	1039.06	0.30	5.76	8.90	15.40	71
2.00	2.01	263.71	1.19	7.73	6.40	15.50	76
3.00	2.98	177.59	1.76	7.73	6.80	14.50	76
3.75	3.67	144.61	2.18	5.54	6.60	9.30	76
4.75	4.87	108.94	2.91	6.40	7.90	10.60	76
5.75	5.86	90.37	3.53	5.76	9.00	9.10	71
6.75	6.79	78.02	4.11	5.76	10.30	9.20	71
8.00	7.91	67.02	4.83	5.76	11.00	8.80	71
9.00	9.10	58.24	5.63	5.69	11.50	8.00	66
11.00	11.00	48.15	6.99	5.76	12.50	5.65	66

Minimum Center Thickness for Surfacing 2.0mm



SHAMIR

shamirlens.com • 877.514.8330