

# Premium FSV Lenses with Glacier PLUS™ AR from Shamir



1.50 Clear / CR-39 Power Range										
SPH	CYL									BASE CURVE
	0.00	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	
0.00										6.00
-0.25										5.50
-0.50										5.00
-0.75										4.75
-1.00										4.50
-1.25										4.50
-1.50										4.25
-1.75										4.00
-2.00										3.75
+0.25										3.00
+0.50										3.00
+0.75										3.25
+1.00										3.50
+1.25										3.75
+1.50										4.00
+1.75										4.25
+2.00										4.50
+2.25										6.25
+2.50										5.00
+2.75										5.25
+3.00										5.50

70mm

65mm

1.53										
SPH	CYL									BASE CURVE
	0.00	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	
0.00										5.00
-0.25										5.00
-0.50										5.00
-0.75										5.00
-1.00										5.00
-1.25										5.00
-1.50										5.00
-1.75										5.00
-2.00										3.00
-2.25										3.00
-2.50										3.00
-2.75										3.00
-3.00										3.00
-3.25										3.00
-3.50										3.00
-3.75										3.00
-4.00										1.00
-4.25										1.00
-4.50										1.00
-4.75										1.00
-5.00										1.00
-5.25										1.00
-5.50										1.00
-5.75										1.00
-6.00										1.00

70mm

65mm

1.59 Polycarbonate Power Range										
SPH	CYL									BASE CURVE
	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	0.00	
0.00										5.00
-0.25										5.00
-0.50										5.00
-0.75										5.00
-1.00										5.00
-1.25										5.00
-1.50										5.00
-1.75										5.00
-2.00										3.00
-2.25										3.00
-2.50										3.00
-2.75										3.00
-3.00										3.00
-3.25										3.00
-3.50										3.00
-3.75										3.00
-4.00										1.00
-4.25										1.00
-4.50										1.00
-4.75										1.00
-5.00										1.00
-5.25										1.00
-5.50										1.00
-5.75										1.00
-6.00										1.00


73mm

70mm\*

65mm

## Are Glacier PLUS™ FSV Lenses suitable for the order?

1. Dot the patients fitting height and PD on the demo lenses of the chosen frame.
2. Place this dot over the chart attached in this packet.
  - If the frame fits the circle, then YES, Shamir Finished Single Vision lenses are a perfect match.
  - If it does not fit, try a Shamir **Freeform** lens like **Autograph III**® Single Vision.
  - If still undecided about which **Freeform** lens to use, call your Shamir Account Executive or **(877) 514-8330** to be connected.



\*available in Blue Zero™ only

1.67 Aspheric Power Range										
SPH	CYL									BASE CURVE
	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	0.00	
0.00										4.00
-0.25										3.75
-0.50										3.50
-0.75										3.25
-1.00										3.00
-1.25										3.25
-1.50										3.00
-1.75										2.75
-2.00										2.50
-2.25										2.25
-2.50										2.00
-2.75										2.25
-3.00										2.00
-3.25										1.75
-3.50										1.50
-3.75										1.75
-4.00										1.50
-4.25										1.25
-4.50										1.00
-4.75										1.25
-5.00										1.00
-5.25										1.25
-5.50										1.00
-5.75										0.75
-6.00										0.50

65mm

75mm

# Finished Single Vision Lenses

Shamir's revolutionary Glacier PLUS™ AR coating is now available on Finished Single Vision lenses!

## **GLACIER PLUS™**

- Shamir Glacier PLUS™ is an Anti-Reflective coating that utilizes advanced technologies to provide a cutting-edge, glare-reducing solution for your patients. It resists scratches, repels water, dust and smudges, and is easier to clean than ever before!

## **ABOUT FINISHED SINGLE VISION**

- Premium quality Shamir branded lenses.
- Faster turn times for patients.
- Available in 1.50 Clear/CR-39, 1.53 Clear, 1.53 Blue Zero™, 1.59 Poly Clear, 1.59 Blue Zero™, & 1.67 Aspheric Clear.

## **COMING SOON**

- **Transitions® GEN8**

(1.59 Gray and Brown)

**NEW  
Product  
Available!**

See Product Availability chart on back for range of lens offerings