


Single Vision HC & UC Spherical



Trivex® NuPolar® Gray & Brown
Certified for Shamir FREEFORM®

Index Nd	1.498
Abbe Value	58
Density gr/cm ³	1.27



NOMINAL BASE CURVE	TRUE BASE (1.530)	BASE RADIUS (MM)	SAG (50MM)	CONCAVE CURVE	CT (MM)	ET (MM)	DIAMETER
2.25	2.26	234.51	1.34	6.25	5.60	10.90	74
4.25	4.24	125.00	2.53	6.25	8.40	11.20	74
6.25	6.25	84.80	3.77	6.25	8.60	8.50	74
8.25	8.26	64.17	5.07	6.25	11.90	8.70	74

Minimum Center Thickness for Surfacing 1.5mm



SHAMIR

shamirlens.com • 877.514.8330