


1.60 Polarized Gray & Brown

Certified for Shamir FREEFORM®



Index Nd	1.664
Abbe Value	31
Density gr/cm ³	1.36

NOMINAL BASE CURVE	TRUE BASE (1.530)	BASE RADIUS (MM)	SAG (50MM)	CONCAVE CURVE	CT (MM)	ET (MM)	DIAMETER	FILM CURVE	FILM LOC @ CTR (MM)
1.75	1.76	301.00	1.040	6.40	9.60	17.50	81	2.49	0.60
2.75	2.73	194.19	1.616	6.40	6.90	13.30	81	2.49	0.9
3.75	3.74	141.69	2.223	6.40	8.00	12.80	81	3.75	0.6
4.75	4.74	111.80	2.831	6.40	8.50	11.70	81	4.45	1.0
5.50	5.50	96.38	3.299	6.40	9.60	11.50	81	5.25	1.0
6.75	6.71	78.98	4.061	6.40	12.10	11.70	81	6.50	1.0
8.50	8.50	62.36	5.231	6.40	14.40	10.40	81	8.23	1.1

Minimum Center Thickness for Surfacing 2.0mm