

Shamir Approved Single Vision Blanks

APPROVED
DESIGNS

AUTOGRAPH INTELLIGENCE™
AUTOGRAPH™ III
AUTOGRAPH™ II+
ATTITUDE III™ - FASHION

ATTITUDE III™ - SPORT
SPECTRUM+™
COMPUTER™
WORKSPACE™

DUO™
AUTOGRAPH™ III - SV
ATTITUDE III™ - SV

	MATERIAL	INDEX	BASE CURVE	COLOR	BLANK SIZE	COATING
HARD RESIN 1.50	Clear	1.5	2, 3.5, 4.25, 6.25, 7.25	Clear	80	HC/UC
	Transitions® Signature™	1.5	0.5, 2, 3.5, 4.25, 5.25, 6.25, 7.25, 8.75, 10.5	Gray/Brown	77	HC/UC
	Transitions® XTRActive™	1.5	0.5, 2.25, 4.25, 5.25, 6.25, 8.25	Gray	76	HC/UC
	Transitions® Vantage™	1.5	0.5, 2.25, 4.25, 5.25, 6.25, 8.25	Gray	76	HC
	Polarized	1.5	2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	Gray/Brown	80	HC
	SunSync®	1.5	2, 4, 5 (GRAY ONLY), 6, 8	Gray/Brown	70/75	HC
	SunSync® Drive XT	1.5	2, 4, 5 (GRAY ONLY), 6, 8	Gray/Brown	70/75	HC
SUPERLIGHT™ 1.60	Clear	1.6	0.5, 2, 3, 4, 5, 6, 7, 8, 9	Clear	75	HC/UC
	Wrap	1.6	4.5, 6, 7, 8, 9, 10.5 (UC ONLY)	Clear	75/80	HC/UC
	Transitions® Signature™	1.6	0.5, 2, 3, 4, 5, 6, 7, 8, 9	Gray/Brown	72	HC/UC
	Transitions® XTRActive™	1.6	0.5, 2, 3, 4, 5, 6, 7, 8, 9	Gray	72	HC/UC
	Polarized	1.6	1.5, 2, 3, 4, 5, 6, 7, 8, 9	Gray/Brown	75	HC
SUPERLIGHT™ 1.67	Clear	1.67	0.5, 2, 3, 4, 5, 6, 7, 8, 9, 11, 12.5	Clear	72	HC/UC
	Blue Zero™	1.67	0.5, 2, 3, 4, 5, 6, 7, 8, 9, 11, 12.5	Blue Zero	75	HC/UC
	Wrap	1.67	4.5, 6, 7, 8, 9, 10 (UC ONLY)	Clear	80	HC/UC
	Transitions® Signature™	1.67	0.5, 2, 3, 4, 5, 6, 7, 8, 9, 11, 12.5	Gray/Brown	72	HC/UC
	Transitions® XTRActive™	1.67	0.5, 2, 3, 4, 5, 6, 7, 8, 9, 11, 12.5	Gray	72	HC/UC
	Polarized	1.67	1.25, 2.75, 3.75, 4.75, 5.25, 6.75, 8	Gray/Brown	81	HC
	SunSync®	1.67	0.5, 2.5, 4, 6.25, 8.5	Gray/Brown	70/75	HC
	SunSync® Drive XT	1.67	0.5, 2.5, 3.25, 4, 4.75, 6.25, 7.75, 8.5	Gray/Brown	70/75	HC
SUPERLIGHT™ 1.74	Clear (HC)	1.74	0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9	Clear	65/72	HC
	Clear (UC)	1.74	0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9	Clear	68/75	UC
	Extra Thick	1.74	0.5, 1, 2	Clear	68/75	UC
	Transitions® Signature™	1.74	0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9	Gray/Brown	65/72	HC/UC
DLC™ TRIVEX™	Clear	1.53	0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9	Clear	75	HC/UC
	Blue Zero™	1.53	0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9	Blue Zero	74	HC/UC
	Transitions® Signature™	1.53	0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9	Gray/Brown	74	HC/UC
	Transitions® Vantage™	1.53	0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9	Gray	74	HC
	Transitions® XTRActive™	1.53	0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9	Gray	74	HC
TRIVEX™	Polarized	1.53	2.25, 4.25, 6.25, 8.25	Gray/Brown	74	HC
	SunSync®	1.53	2, 4, 6, 8	Gray/Brown	70/75	HC
	SunSync® Drive XT	1.53	2, 4, 6, 8	Gray/Brown	70/75	HC
POLYCARBONATE	Clear	1.59	0.75, 1.5, 3.25, 4.25, 5.5, 6.5, 7.5, 8.5, 9.75	Clear	76/81	HC
	Blue Zero™	1.59	0.5, 1.5, 2.0, 3.5, 4.25, 5.25, 6, 7.5, 8.5	Blue Zero	75	HC
	Transitions® Signature™	1.59	0.75, 1.5, 3.25, 4.25, 5.5, 6.5, 7.5, 8.5, 9.75	Gray/Brown	74/79	HC
	Transitions® Vantage™	1.59	0.75, 1.5, 3.25, 4.25, 5.5, 6.5, 7.5, 8.5, 9.75	Gray	74/79	HC
	Transitions® XTRActive™	1.59	0.75, 1.5, 3.25, 4.25, 5.5, 6.5, 7.5, 8.5, 9.75	Gray	74/79	HC
	Polarized	1.59	1.5, 2.25, 3.25, 4.25, 5.5, 6.5, 7.5, 8.5	Gray/Brown	76	HC
	Polarized (Wrap)	1.59	6.25, 7.25, 8.25	Gray/Brown	82	HC
	Transitions® Drivewear®	1.59	2.25, 4.25, 6.25, 8.25		75	HC
	SunSync®	1.59	0.5, 2, 4, 6, 8	Gray/Brown	70/75	HC
	SunSync® Drive XT	1.59	0.5, 2, 4, 6, 8	Gray/Brown	70/75	HC
NXT	NXT	1.53	6, 8	Green/Copper*	80	HC
	Polarized	1.53	8	Gray/Brown*	80	HC

OTHER AUTHORIZED MANUFACTURERS	INDEX	MATERIAL	MANUFACTURER	INDEX	MATERIAL	MANUFACTURER
	1.5	1.5	Polarized G-15	KBCO	1.59	BluTech® Ultra Indoor
Transitions® Drivewear®			Younger Optics	BluTech® Ultra Outdoor Polarized		BluTech
Transitions® Signature™ Graphite Green			Essilor	Polarized G-15		KBCO
Transitions® XTRActive™			Essilor	Transitions® Signature™ Graphite Green		Essilor
Xperio Polarized			Essilor	Xperio Polarized		Essilor
1.53	1.53	Trivex™ Polarized G-15	KBCO & X-Cel	1.6	Polarized G-15	KBCO
		Xperio Polarized	Essilor		Xperio Polarized	Essilor
1.56	1.56	BluTech® Indoor	BluTech	1.67	Polarized G-15	KBCO
		BluTech® Outdoor Polarized	BluTech		Transitions® Signature™ Graphite Green	Essilor
					Xperio Polarized	Essilor

* WHILE SUPPLIES LAST (LIMITED QUANTITIES IN STOCK)
RED BASE CURVE COINCIDES WITH RED BLANK SIZE



Base Curves, Power Range, etc. can be obtained from your lens supplier.