


1.67 Transitions® Gen 8™ Gray & Brown Certified for Shamir FREEFORM®

Index Nd	1.665
Abbe Value	32
Density gr/cm ³	1.35



NOMINAL BASE CURVE	TRUE BASE (1.530)	BASE RADIUS (MM)	SAG (50MM)	CONCAVE CURVE	CT (MM)	ET (MM)	DIAMETER
0.50	0.50	1060.00	0.29	5.76	8.90	15.40	71
2.00	1.99	265.82	1.18	7.73	6.40	15.50	76
3.00	2.97	178.54	1.75	7.73	6.80	14.50	76
3.75	3.65	145.20	2.17	5.54	6.60	9.30	76
4.75	4.85	109.30	2.90	6.40	7.90	10.60	76
5.75	5.85	90.62	3.51	5.76	9.00	9.10	71
6.75	6.78	78.20	4.10	5.76	10.30	9.20	71
8.00	7.89	67.16	4.82	5.76	11.00	8.80	71
9.00	9.08	58.34	5.62	5.69	11.50	8.00	66
11.00	11.00	48.18	6.99	5.76	12.50	5.65	66

Minimum Center Thickness for Surfacing 2.0mm