

**WORK**

Shamir

**Autograph II - Office™**

**PRESCRIPTOR ORDER CODE:**







AUTOGRAPH II OFFICE








**PRODUCT FEATURES:**







- FREEFORM PAL
- FREEFORM OCCUPATION
- FREEFORM SV
- CONVENTIONAL PAL
- CONVENTIONAL OCCUPATION
- CONVENTIONAL SV
- EYEPOINT TECHNOLOGY®
- AS-WORN TECHNOLOGY®







**NOTE:**

- Recommended min fitting height is 16mm
- Prism to 10.00 Diopters
- Maximum dynamic reduction -2.25

CR-39	TREATMENT		COATING
	☺ Clear		
TOTAL NEAR POWER RANGE	-6.00 to +6.00		
PRESCRIBED ADD POWER RANGE	0.75 to 3.50		
CYLINDER	to -6.00		
OTHER			
INVISIBLE MARKINGS	 -175  15  FFO	 BZP  FFO	

POLYCARBONATE	TREATMENT		COATING
	☺ Clear	☺ Blue Zero™	
TOTAL NEAR POWER RANGE	-6.00 to +8.00	-6.00 to +8.00	
PRESCRIBED ADD POWER RANGE	0.75 to 3.50	0.75 to 3.50	
CYLINDER	to -7.00	to -7.00	
OTHER			
INVISIBLE MARKINGS	 -175  PC  FFO	 -175  BZP  FFO	

SuperLite™ 1.60	TREATMENT		COATING
	☺ Clear		
TOTAL NEAR POWER RANGE	-6.00 to +8.00		
PRESCRIBED ADD POWER RANGE	0.75 to 3.50		
CYLINDER	to -7.00		
OTHER			
INVISIBLE MARKINGS	 -175  16  FFO	 BZP  FFO	

 Shamir Glacier PLUS™ UV
  Crizal AR Coatings
  Hoya AR Coatings
  iCoat AR Coatings
  Zeiss AR Coatings
  Unity AR Coatings

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

☺ - Denotes that material is available in a Shamir Freeform Blank™

**EFFECTIVE DATE: FEBRUARY 01, 2018**

