


# Shamir Autograph II - Attitude™

EFFECTIVE DATE: NOVEMBER 01, 2020



**NOTES**

- Recommended min. fitting heights are 15mm and 18mm
- Prism to 10.00 Diopters

**PRESCRIBTOR ORDER CODE**

AUTOGRAPH II ATTITUDE 15  
AUTOGRAPH II ATTITUDE 18

**PRODUCT FEATURES**

- SHAMIR AI ENGINE™
  - AS-WORN QUADRO™
  - AS-WORN TECHNOLOGY™
  - CLOSE-UP™
  - CONTINUOUS DESIGN TECHNOLOGY™
  - EYE-POINT TECHNOLOGY®
  - EYE-POINT TECHNOLOGY AI™
  - EYE-POINT TECHNOLOGY III®
  - FREEFORM PAL
  - FREEFORM SV
  - INTELICORRIDOR™
  - NATURAL POSTURE™
- PATENT PENDING

		MATERIAL / TREATMENT													
		HARD RESIN 1.50		TRIVEX™		1.56		POLYCARBONATE 1.59		1.60		1.67		1.74	
		POWER RANGE	-5.00 to +4.00	POWER RANGE	-5.00 to +4.00	POWER RANGE	-5.50 to +6.00	POWER RANGE	-6.50 to +6.00	POWER RANGE	-6.50 to +6.00	POWER RANGE	-6.50 to +5.00	POWER RANGE	
		ADD POWER	+0.75 to +3.50	ADD POWER	+0.75 to +3.50	ADD POWER	+0.75 to +3.50	ADD POWER	+0.75 to +3.50	ADD POWER	+0.75 to +3.50	ADD POWER	+0.75 to +3.50	ADD POWER	
		CYLINDER	to -6.00	CYLINDER	to -6.50	CYLINDER	to -6.75	CYLINDER	to -7.00	CYLINDER	to -7.00	CYLINDER	to -8.00	CYLINDER	
CLEAR	15	TX	56	PC	16	67									
	TR	TTX		PT	T6	T67									
Transitions®	GRAY/BROWN	GRAY/BROWN		GRAY/BROWN	GRAY/BROWN	GRAY/BROWN		GRAY/BROWN	GRAY/BROWN	GRAY/BROWN	GRAY/BROWN	GRAY/BROWN	GRAY/BROWN		
XTRACTIVE®	X5	XX		PX	X6	X67									
	GRAY	GRAY		GRAY	GRAY	GRAY		GRAY	GRAY	GRAY	GRAY	GRAY	GRAY		
sunsync	S50	S53		SPY		S67									
VANTAGE®	V5	VX		PV											
	GRAY	GRAY		GRAY											
POLARIZED	PZ	PTX	56P	PPZ	PZ6	PZ67									
	GRAY/BROWN G-15	GRAY/BROWN G-15	BLUTECH OUTDOOR	GRAY/BROWN G-15	GRAY/BROWN G-15	GRAY/BROWN G-15		GRAY/BROWN G-15	GRAY/BROWN G-15	GRAY/BROWN G-15	GRAY/BROWN G-15	GRAY/BROWN G-15	GRAY/BROWN G-15		
DRIVEWEAR®	DW			PD											
BLUE ZERO™		BZX		BZP		BZ67									
NXT		NXTP													
		POWER RANGE	-5.00 to +4.00												
		GREEN/COPPER													

- LEGEND**
- Shamir Glacier PLUS™
  - Crizal AR Coatings
  - Hoya AR Coatings
  - iCoat AR Coatings
  - Unity AR Coatings
  - Zeiss AR Coatings

\* Photochromic solutions include Transitions®, XTRACTive®, SunSync®, SunSync® Drive XT (GRAY ONLY), SunSync® Elite, & SunSync® Elite XT.  
 \*\* PATENT PENDING

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

**⚠ WARNING:** Polycarbonate lenses can expose you to chemicals including Bisphenol A (BPA) and 1.74 index lenses can expose you to Methimazole, which are known to the State of California to cause birth defects or other reproductive harm. For more information go to: [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

SHI-A2-ATT-REF-111120