

Shamir Autograph II - Office™

EFFECTIVE DATE: NOVEMBER 01, 2020

NOTES


- Recommended min fitting height is 16mm
- Prism to 10.00 Diopters
- Maximum dynamic reduction -2.25

**PRESCRIBTOR
ORDER CODE**

AUTOGRAPH II OFFICE

PRODUCT FEATURES

- SHAMIR AI ENGINE™
- AS-WORN QUADRO™
- AS-WORN TECHNOLOGY™
- CLOSE-UP™
- CONTINUOUS DESIGN TECHNOLOGY™
- EYE-POINT TECHNOLOGY™
- EYE-POINT TECHNOLOGY AI™
- EYE-POINT TECHNOLOGY III™
- FREEFORM OCCUPATION
- FREEFORM SV
- INTELLICORRIDOR™
- NATURAL POSTURE™



	MATERIAL / TREATMENT													
	HARD RESIN 1.50		TRIVEX™		1.56		POLYCARBONATE 1.59		1.60		1.67		1.74	
	POWER RANGE	-6.00 to +6.00	POWER RANGE		POWER RANGE		POWER RANGE	-6.00 to +8.00	POWER RANGE	-6.00 to +8.00	POWER RANGE		POWER RANGE	
	ADD POWER	+0.75 to +3.50	ADD POWER		ADD POWER		ADD POWER	+0.75 to +3.50	ADD POWER	+0.75 to +3.50	ADD POWER		ADD POWER	
	CYLINDER	to -6.00	CYLINDER		CYLINDER		CYLINDER	to -7.00	CYLINDER	to -7.00	CYLINDER		CYLINDER	
CLEAR	15						PC		16					
Transitions®														
XTRACTIVE®														
sunsync														
VANTAGE®														
POLARIZED														
DRIVEWEAR®														
BLUE ZERO™							BZP							
NXT														

- LEGEND**
- Shamir Glacier PLUS™
 - Crizal AR Coatings
 - Hoya AR Coatings
 - iCoat AR Coatings
 - Unity AR Coatings
 - Zeiss AR Coatings

* Photochromic solutions include Transitions®, XTRActive®, SunSync®, SunSync® Drive XT (GRAY ONLY), SunSync® Elite, & SunSync® Elite XT.
 ** PATENT PENDING

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

⚠ **WARNING:** Polycarbonate lenses can expose you to chemicals including Bisphenol A (BPA) and 1.74 index lenses can expose you to Methimazole, which are known to the State of California to cause birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov.