

Shamir Autograph II - Single Vision Attitude™

EFFECTIVE DATE: NOVEMBER 01, 2020

NOTES


• Prism to 10.00 Diopters

PRESCRIBTOR ORDER CODE

AUTOGRAPH II SV ATTITUDE

PRODUCT FEATURES

- SHAMIR AI ENGINE™
- AS-WORN QUADRO™
- AS-WORN TECHNOLOGY™
- CLOSE-UP™
- CONTINUOUS DESIGN TECHNOLOGY™
- EYE-POINT TECHNOLOGY®
- EYE-POINT TECHNOLOGY AI™
- EYE-POINT TECHNOLOGY III®
- FREEFORM SV
- FREEFORM PAL
- INTELICORRIDOR™
- NATURAL POSTURE™



	MATERIAL / TREATMENT													
	HARD RESIN 1.50		TRIVEX™		1.56		POLYCARBONATE 1.59		1.60		1.67		1.74	
	POWER RANGE	ADD POWER	POWER RANGE	ADD POWER	POWER RANGE	ADD POWER	POWER RANGE	ADD POWER	POWER RANGE	ADD POWER	POWER RANGE	ADD POWER	POWER RANGE	ADD POWER
	-6.00 to +6.00		-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-20.00 to +15.50	-20.00 to +15.50	-20.00 to +15.50
	CYLINDER		CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER
	to -6.00		to -6.50	to -6.75	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -8.00	to -8.00	to -9.00	to -9.00	to -9.00
CLEAR	15		TX		56		PC		16		67		75	
Transitions®	TR		TTX				PT		T6		T67		T74	
XTRACTIVE®	X5		XX				PX		X6		X67			
sunsync	S50		S53				SPY				S67			
VANTAGE®	V5		VX				PV							
POLARIZED	PZ		PTX		56P		PPZ		PZ6		PZ67			
DRIVEWEAR®	DW						PD							
BLUE ZERO™			BZX				BZP				BZ67			
NXT			NXTP											
			POWER RANGE											
			-6.00 to +6.00											
			GREEN/COPPER											

- LEGEND**
- Shamir Glacier PLUS™
 - Crizal AR Coatings
 - Hoya AR Coatings
 - iCoat AR Coatings
 - Unity AR Coatings
 - Zeiss AR Coatings

* Photochromic solutions include Transitions®, XTRActive®, SunSync®, SunSync® Drive XT (GRAY ONLY), SunSync® Elite, & SunSync® Elite XT.
 ** PATENT PENDING

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

⚠ **WARNING:** Polycarbonate lenses can expose you to chemicals including Bisphenol A (BPA) and 1,74 index lenses can expose you to Methimazole, which are known to the State of California to cause birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov.

SHM-A23V-ATT-REF-11120