

EFFECTIVE DATE: NOVEMBER 01, 2020

NOTES


- Recommended min. fitting height is 15mm
- Prism to 10.00 Diopters

PRESCRIBTOR ORDER CODE

DUO

PRODUCT FEATURES

- SHAMIR AI ENGINE™
- AS-WORN QUADRO™
- AS-WORN TECHNOLOGY™
- CLOSE-UP™
- CONTINUOUS DESIGN TECHNOLOGY™
- EYE-POINT TECHNOLOGY®
- EYE-POINT TECHNOLOGY AI™
- EYE-POINT TECHNOLOGY III®
- FREEFORM PAL
- FREEFORM BIFOCAL
- INTELLICORRIDOR™
- NATURAL POSTURE™



		MATERIAL /TREATMENT						
		HARD RESIN 1.50	TRIVEX™	1.56	POLYCARBONATE 1.59	1.60	1.67	1.74
		POWER RANGE -12.25 to +6.00	POWER RANGE -13.25 to +6.00	POWER RANGE -14.00 to +6.00	POWER RANGE -14.75 to +6.00	POWER RANGE -15.00 to +6.00	POWER RANGE -16.75 to +8.00	POWER RANGE -18.75 to +12.00
		ADD POWER +0.75 to +4.00						
		CYLINDER to -6.00	CYLINDER to -6.50	CYLINDER to -6.75	CYLINDER to -7.00	CYLINDER to -7.00	CYLINDER to -8.00	CYLINDER to -9.00
CLEAR	15	TX	56	PC	16	67	74	
Transitions®	TR	TTX		PT	T6	T67	T74	
XTRACTIVE®	X5	XX		PX	X6	X67		
sunsync	S50	S53		SPY		S67		
VANTAGE®	V5	VX		PV				
POLARIZED	PZ	PTX	56P	PPZ	PZ6	PZ67		
DRIVEWEAR®	DW			PD				
BLUE ZERO™		BZX		BZP		BZ67		
NXT		NXTP						
		POWER RANGE 1.75 to +6.00						
		GREEN/COPPER						

- LEGEND**
- Shamir Glacier PLUS™
 - Crizal AR Coatings
 - Hoya AR Coatings
 - iCoat AR Coatings
 - Unity AR Coatings
 - Zeiss AR Coatings

* Photochromic solutions include Transitions®, XTRActive®, SunSync®, SunSync® Drive XT (GRAY ONLY), SunSync® Elite, & SunSync® Elite XT.
 ** PATENT PENDING

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

⚠ **WARNING:** Polycarbonate lenses can expose you to chemicals including Bisphenol A (BPA) and 1.74 index lenses can expose you to Methimazole, which are known to the State of California to cause birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov.

SHM-DUO-REF-11120