

NOTES

- Recommended min fitting heights are 14mm, 16mm, and 18mm
- Prism to 10.00 Diopters


PRESCRIBTOR ORDER CODE

- SPECTRUM 14
- SPECTRUM 16
- SPECTRUM 18

PRODUCT FEATURES

- SHAMIR AI ENGINE™
- AS-WORN QUADRO™
- AS-WORN TECHNOLOGY™
- CLOSE-UP™
- CONTINUOUS DESIGN TECHNOLOGY™
- EYE-POINT TECHNOLOGY®
- EYE-POINT TECHNOLOGY AI™
- EYE-POINT TECHNOLOGY III®
- FREEFORM PAL
- FREEFORM SV
- INTELICORRIDOR™
- NATURAL POSTURE™

PATENT PENDING



		MATERIAL / TREATMENT						
		HARD RESIN 1.50	TRIVEX™	1.56	POLYCARBONATE 1.59	1.60	1.67	1.74
		POWER RANGE	POWER RANGE	POWER RANGE	POWER RANGE	POWER RANGE	POWER RANGE	POWER RANGE
		-12.25 to +6.00	-13.25 to +6.00	-14.00 to +6.00	-14.75 to +6.00	-15.00 to +6.00	-16.75 to +8.00	-18.75 to +14.74
		ADD POWER	ADD POWER	ADD POWER	ADD POWER	ADD POWER	ADD POWER	ADD POWER
		+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50
		CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER	CYLINDER
		to -6.00	to -6.50	to -6.75	to -7.00	to -7.00	to -8.00	to -9.00
CLEAR	15	TX	56	PC	16	67	74	
	TR	TTX		PT	T6	T67	T74	
Transitions®	X5	XX		PX	X6	X67		
	S50	S53		SPY		S67		
XTRACTIVE®	V5	VX		PV				
	PZ	PTX	56P	PPZ	PZ6	PZ67		
sunsync	DW			PD				
		BZX		BZP		BZ67		
VANTAGE®								
POLARIZED								
DRIVEWEAR®								
BLUE ZERO™								
NXT								

- LEGEND**
- Shamir Glacier PLUS™
 - Crizal AR Coatings
 - Hoya AR Coatings
 - iCoat AR Coatings
 - Unity AR Coatings
 - Zeiss AR Coatings

* Photochromic solutions include Transitions®, XTRACTIVE®, SunSync®, SunSync® Drive XT (GRAY ONLY), SunSync® Elite, & SunSync® Elite XT.
 ** PATENT PENDING
 All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.
 ⚠ **WARNING:** Polycarbonate lenses can expose you to chemicals including Bisphenol A (BPA) and 1.74 index lenses can expose you to Methimazole, which are known to the State of California to cause birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov.